



# Levo Plus Motorized Pipette Filler

## Levo Plus





Levo Plus c/w AC adapter, spare 0.45um filter and wall stand

### Features:

- ❑ A full volume range of 1-100ml
- ❑ Easy to one-handed operation
- ❑ Lightweight, ergonomic design allows longer, fatigue-free pipetting
- ❑ Suitable for all types of pipettes 1-100ml
- ❑ Lithium-ion battery provides long runtime charge and no battery memory problems
- ❑ Low battery indication
- ❑ Powerful motor fills a 25ml pipette in <5seconds
- ❑ Supplied with integral but replaceable 0.45um or 0.2um filter
- ❑ Can be used while recharging
- ❑ **3 YEAR WARRANTY**

### Ordering information

Cat. No.	Description
740200019999	Levo Plus c/w 110V AC Adapter & spare 0.45um, filter & wall stand

### Accessories

Cat. No.	Description
17000103	0.45um Hydrophobic Filter, 1ea/PK
17000144	0.2um Hydrophobic Filter, 1ea/PK

### Main Data Description

Aspirate Speeds	8 coarse settings plus fine speed control with triggers
Dispense Speeds	8 & Gravity Dispense
Battery	Lithium-Ion
Battery Service Life	More than 8 Hours
Charging Time	Approx.2 Hours
Pipette Types	Glass or Plastic Pipettes 1-100ml
Filter	0.2µm or 0.45µm Hydrophobic

**DRAGONLAB**

500 Four Rod Road  
Suite 122  
Berlin, CT 06037

Tel: (860) 828-5289  
Fax: (860) 828-5389  
info@dragonlabusa.com  
www.dragon-lab.com